



Digital Neurotherapeutics

**Neurotechnology for everybody, anytime
by anyone**

Javier Mínguez

Scientific Director of Bitbrain

Professor Universidad de
Zaragoza

@bitbrain_en | www.bitbrain.com



First steps in neurotechnology (2003-2012)



Second demonstration of wheelchair control only with brain thoughts
One of the most cited papers on BCI for robot control
Best Science & Technology video of the month (April 2009, Global)



Universidad
Zaragoza

- Pioneering brain-computer interfaces (BCI)
- 150+ scientific papers
- Large scientific impact in Robotics and BCI
- Very large traction in EU projects.

Research featured in:

nature

WIRED


BBC

Science

NewScientist

The
Economist


First steps in neurotechnology (2003-2012)




First demonstration of a Robot that learns its operation from brain signals (2015)



Universidad
Zaragoza




First steps in neurotechnology (2003-2012)



First demonstration that natural motor intention can be reliably captured in stroke patients from healthy motor cortex




Universidad
Zaragoza



Mission

Bringing neurotechnology closer to society



I+D requerida

- EEG systems and Sw
- EEG SP and IA
- EEG biomarkers
- BCI treatments

Others

- Regulatoria
- Seguridad
- Etica
- Stakeholders



2003



Universidad
Zaragoza

2012



2023

Real world neurotech

Research company

60 people with organic growth.

R&D and Technology

- A team of more than 20 engineers and 15 R&D researchers (and growing!).
- Significant investment in R&D (+€6 million invested to date)

Commercialization

Distributors in BRICS countries (Brazil, Russia, India, China and South Africa).

Clients in more than 60 countries (Europe, US, Asia, Latam & Africa)




What is Digital Neurotherapeutics?

Addressing the Continuum of Neurology and Mental Health

Digital neurotherapeutics

Prevent, treat, and manage with software-based solutions


Overcome the limitations of time, place, and a person's reach in delivering healthcare.



EEG recording in Digital Neurotherapeutics


Wearable garments that measure brain activity (medical device)

Equipment usability



A 12-year-old child using the equipment autonomously

Technical evaluation wrt medical-grade EEG



- [Paper](#): Evaluation in Spontaneous and Evoked EEG Time-Frequency.
- Ongoing research for a medical device, such as self-administered partial EEG and PSG (with PPG) – for use by laypersons at home.

AI for sleep scoring


Automatic sleep labeling based on EEG (medical device)

Estado del arte

- Reliability is 80-85% among expert evaluators when manual scoring is performed [4].
- 83% with full PSG and deep learning AI [4].
- 82% with only two frontal sensors and deep learning AI [5].
- 83% (Bitbrain) replicated [5] with N=197 (Base SleepEDF) [6].

Objetivo y estado actual

- Develop an AI that achieves >80% accuracy with data from the band (frontal sensors).
- Currently, with a sample size of N=100, it achieves 87% with PSG and 84% with the band.



[4] Inter-rater reliability when performing manual scoring is 80%, meaning that two well-trained scorers will, on average, agree on 80% of the labels they assign to each epoch.

[5] A. Koushik et al., Real-Time Sleep Staging using Deep Learning on a Smartphone for a Wearable EEG, Machine Learning for Health (ML4H) Workshop at NeurIPS 2018

[6] M. Esparza-laizzo et al. Automatic sleep scoring for real-time monitoring and stimulation in individuals with and without sleep apnea. BioRxiv. 2024

Cognitive deterioration in DTx

HOGAR: MCI/ALZ Monitorization


Estudio HOGAR -- 2024

- **Objective:** Monitoring the continuum of cognitive decline through self-administered EEG and PSG at the patient's home.
- **Population:** 500 participants across the Alzheimer's continuum Healthy, Subjective memory complaints, Mild Cognitive Impairment (MCI), Mild Alzheimer's Disease.



Status



- **Status:** 103/200 (80% success rate)
- Home Instrument Results:
 - Self-administered EEG: 100% (2 recordings per subject ~15 minutes)
 - Self-administered PSG: 100% (2 recordings per subject: ~8 hours each)



Treatment of Cognitive Deterioration in DTx (I)


Neurological basis

Greater desynchronization of individualized high-alpha rhythms in the parietal and parieto-occipital regions is associated with higher general cognitive ability, particularly in sustained attention and working memory.




Objective:

Use neurofeedback to induce a change in desynchronization during task-related EEG.




Treatment of Cognitive Deterioration in DTx (II)

Neuromodulation through Neurofeedback (active treatment)



Treatment of Cognitive Deterioration in DTx (III)

Evidence and current studies



References

- Escolano, C., Aguilar, M., & Mínguez, J. (2011). EEG-based upper alpha neurofeedback training improves working memory performance. In International Conference of the IEEE Engineering in Medicine and Biology Society (EMBC) (p. 2327–2330). Boston (USA).
- Escolano, C., Navarro-Gil, M., García-Campayo, J., Congedo, M., De Ridder, D., & Mínguez, J. (2014a). A controlled study on the cognitive effect of alpha neurofeedback training in patients with major depressive disorder. *Frontiers in Behavioral Neuroscience*, 8(296).
- Escolano, C., Navarro-Gil, M., García-Campayo, J., Congedo, M., & Mínguez, J. (2014b). The effects of individual upper alpha neurofeedback in ADHD: an open-label pilot study. *Applied psychophysiology and biofeedback*, 39, 193–202.
- Navarro-Gil, M., Escolano, C., Montero-Marín, J., Mínguez, J., Shonin, E., & García-Campayo, J. (2018). Efficacy of neurofeedback on the increase of mindfulness-related capacities in healthy individuals: a controlled trial. *Mindfulness*, 9, 303–311.



- Preclínica ESPERANZA (2023): doble ciego controlado randomizado N=20 de 40 a 60 años. 10 sesiones de 1h en 2 semanas.
- Feasibility ESPERANZA (2023-2024): doble ciego controlado randomizado con N=30 DCLa de 30 sesiones de 1h en 1,5 meses.
- Pivotal ESPERANZA (2023-2024): doble ciego controlado randomizado con N=40 DCLa de 30 sesiones de 1h en 1,5 meses.

Treatment of Cognitive Deterioration in DTx (I)

Neurological basis



Increasing the number and magnitude of slow waves during deep sleep improves declarative memory.

During the depolarized phase (positive surface) of slow waves in the EEG during Non-REM sleep, well-calibrated auditory stimulation can produce slow oscillations without causing awakening.



Objetivo :

Utilizar la estimulación
auditiva durante el
sueño para inducir
trenes de ondas
lentas adicionales




[8] Ngo et al. Ngo HV, Martinetz T, Born J, Mölle M. Auditory closed-loop stimulation of the sleep slow oscillation enhances memory. *Neuron*. 2013

[9] J. G. Klinzing et al., Mechanisms of systems memory consolidation during sleep. *Nature neuroscience*. 2019

Treatment of Cognitive Deterioration in DTx (II)

Neuromodulation through auditory stimulation during sleep

Neuromodulation by auditory stimulation during sleep for cognitive enhancement




IA EEG filtering
(sleep staging)



EEG signal processing and auditory stimulation

OFF



Treatment of Cognitive Deterioration in DTx (II)


Current Evidence and Project Status

Auditory stimulation (cognitive outcome)

Hedges' g 95 % CL

Ngo et al., 2013
Ngo et al., 2015
Henin et al., 2019 Exp.1
Henin et al., 2019 Exp.2
Ong et al., 2016
Schneider et al., 2020
Papalambros et al., 2017
Leminen et al., 2017

combined effect size



[10] Wunderlin, M. et al. Modulating overnight memory consolidation by acoustic stimulation during slow-wave sleep: a systematic review and meta-analysis. *Sleep*, 44(7). 2021.

- Preclínica NANA (2023): estudio de investigación N=26 personas de 25 a 81 años. 1-3 noches a la semana.
- Feasibility NANA (2023-2024): doble ciego controlado randomizado con N=10 DCLa de 30 sesiones de 1h en 1,5 meses.
- Pivotal NANA (2023-2024): doble ciego controlado randomizado con N=37 DCLa de 30 sesiones de 1h en 1,5 meses.

Participa:    

Financia:   

Aprobado: CEICA (Ref: C.I. PI 23/186 y C.I. EC 22/011 | AEMPS (Ref: 1040/22/EC-R)

Conclusion

Digital neurotherapeutics

Take home messages

- Prevent, manage, and treat with software-based solutions
- Overcome the limitations of time, location, and accessibility in delivering healthcare.
- Development of low-density EEG and PSG sensors for use by laypersons, allowing these records to be taken outside of healthcare facilities.
- Central role of artificial intelligence in automatically analyzing data and personalizing interventions.
- Approach enables very intensive treatments both during wakefulness and sleep, which are impractical to implement in healthcare centers.



The team!






Advanced
neurotechnology


Participantes/Voluntarios

neurociencia@bitbrain.es


www.bitbrain.com

Find us:

 @Bitbrain_EN | @Bitbrain_ES

 @BitbrainEN | @BitbrainES

 @Bitbrain_EN

 Bitbrain

 Bitbrain